


# 2018 RME Crossword Puzzle



**Across**

- 3. Equipment that is used repeatedly is considered \_\_\_\_\_ medical equipment and needs to be cleaned between each new patient use.
- 5. All reusable instruments must be opened/unhinged and sprayed with \_\_\_\_\_ cleaner before sending them to sterile processing.
- 8. You must use Clorox \_\_\_\_\_ to clean the 12-lead EKG machines
- 10. Hydrogen peroxide has a \_\_\_\_\_ minute wet time; unless TB was suspected, then it is five minutes.
- 11. Wipe equipment down from the \_\_\_\_\_ areas to the dirtiest areas.

**Down**

- 1. Stainless steel instruments that are commonly shiny and are stamped with sole use, Pakistan or 3M, should be disposed of in the \_\_\_\_\_ container.

- 2. Commode chairs that are used by patients on contact isolation must \_\_\_\_\_ in the patient room for EMS staff to clean.
- 4. Remove any visibly \_\_\_\_\_ with a disposable cloth or gauze pad before cleaning equipment or spraying down instruments.
- 6. If an instrument pack has been damaged or if sterility has been compromised it should \_\_\_\_\_ be used.
- 7. External pacemakers (epicardial/transvenous) should ONLY be cleaned with \_\_\_\_\_.
- 9. When using sterilized instrument packs you must check that it was properly sterilized by looking for a \_\_\_\_\_ change on the chemical indicator tape BEFORE opening the package for use.