


Chest Tube Assessment



Across

6. The orange suction bellow should be expanded to the _____ mark

7. the patient pressure float ball will fall with _____

8. When assessing for an air leak, you assess for _____ to be present in the air leak chamber

9. The drainage system should always be _____ the level of the patient's chest

11. Chest tube tubing should or should not have dependent loops

12. Tidaling is observed in the _____ seal chamber

Down

1. Tidaling will _____ as the pneumo decreases in size

2. The wall suction should be set at _____ mmHg

3. Changing suction pressure is accomplished by adjusting the rotary of the Dry Suction _____ on the upper left hand side of the Oasis

4. The orange accordian is called the Suction Monitor _____

5. The patient pressure float ball will rise with _____

10. Our policy is to assess the chest tube managment system every _____ hours