

Name: _____ Date: _____

TOOLS OF THE TRADE

1. WHAT COLOR TUBE DO YOU DRAW A PROTHROMBIN TIME IN ? A. 22
2. WHAT SIZE NEED IS USED TO DRAW A DONOR UNIT OF BLOOD? B. BLUE
3. WHAT HOLDS A DOUBLE POINTED NEEDLE? C. SST
4. WHAT IS SHARP ON BOTH ENDS? D. DOUBLE POINTED NEEDLE
5. USED TO SLOW THE RETURN OF VENOUS BLOOD TO THE HEART E. TENDERLET
6. WHAT YOU DISCARD NEEDLES IN? F. TAN
7. WHERE YOU SIT TO HAVE BLOOD DRAWN. G. SODIUM
8. HAS WINGS. H. NEEDLE HOLDER
9. WHAT HAS SODIUM POLYANETHOL SULFONATE (SPS} ADDITIVE IN IT? I. 18
10. USED TO DO SKIN PUNCTURES ON BABIES. J. PINK
11. DO NOT USE A GREEN TUBE WITH SODIUM HEPARIN FOR THIS TEST. K. BLOOD CULTURE BOTTLE
12. WHAT COLOR TUBE HAS EDTA AND USED FOR CBC'S? L. BUTTERFLY
13. WHAT TUBE IS USED TO SEPARATE CELLS AND SERUM WHEN CENTRIFUGED? M. TRANSFER UNIT
14. WHAT COLOR TUBE IS USED TO DO A BLOOD TYPE? N. SHARPS CONTAINER
15. WHAT COLOR TUBE IS USED FOR LEAD TESTING? O. PHLEBOTOMY CHAIR
16. WHAT IS THE NEXT SMALLER NEEDLE TO A 21 GAUGE? P. LAVENDER
17. WHAT IS USED TO TRANSFER BLOOD FROM A SYRINGE TO A TUBE? Q. TOURNIQUET